

## RESECTOSCOPY SHEATH

13Fr.; 1 turnable stopcock

pediatric



**SKU:** 630-510-13

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Pediatric Resectoscopy</a>
<b>TYPE</b>	<a href="#">Sheaths</a>