

RETRIEVAL BASKET, STAINLESS

hexagonal; 20 X 60mm; 4 wires

2.5mm.; 210cm; Lithotripsy



SKU: PK-1SHPK2520210AK

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)