

RETRIEVAL BASKET, STAINLESS

spiral; 20 x 40mm; wires

2.5mm; 210cm, Lithotripsy



SKU: PK-1SSPK2520210VAK

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)