

## RETRIEVAL BASKET, STAINLESS

hexagonal, 20 x 40mm, 4 wires, f. 0,035"

2.5mm, 230cm



**SKU:** PK-1FHPK2520230K

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Stonebaskets for Retrieval & Lithotripsy](#)