

RETRIEVAL BASKET, STAINLESS

hexagonal, 30 x 60mm, 4 wires, f. 0,035"

2.5mm, 230cm



SKU: PK-1FHPK2530230K

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)