

RETRIEVAL BASKET, STAINLESS

hexagonal, 10 x 20mm, 4 wires

2.3mm; 230cm



SKU: PK-1BHPK2310230K

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)