

RING APPLICATOR FOR 7.00MM

tube sterilization

Ø 6.5mm; 300mm; trocare

SKU: 587-230-30
Categories: [Endoscopy](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Laparoscopy Instruments
TYPE	Ring Aplicators

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Laparoscopy Instruments
TYPE	Ring Aplicators