

RING STRIPPER, ROUND MODEL

Ø 4.0mm

50.0cm (20")

SKU: 080-420-04

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|------------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Strippers |
| KLS Martin | 16-251-04 |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|------------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Strippers |
| KLS Martin | 16-251-04 |