

ROCHESTER-PEAN HAEMOSTATIC FORCEPS

straight; serrated

30.0cm (12")



SKU: 040-260-30
Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH474
KLS Martin	13-410-30
AUTHORS	PEAN , ROCHESTER
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps
Aesculap	BH474
KLS Martin	13-410-30
AUTHORS	PEAN , ROCHESTER

SPECIFICATIONS

[serrated](#), [straight](#)