

## ROGERS HYSTERECTOMY FORCEPS

straight; serrated

25.5cm (10")



**SKU:** 040-580-24

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ550
<b>KLS Martin</b>	13-637-10
<b>AUTHORS</b>	<a href="#">ROGERS</a>
<b>Medicon</b>	16.35.50
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Hysterectomy Forceps</a>
<b>Aesculap</b>	BJ550
<b>KLS Martin</b>	13-637-10

## AUTHORS

[ROGERS](#)

## Medicon

16.35.50

## SPECIFICATIONS

[serrated](#), [straight](#)