

## RUMEL-BELMONT TOURNIQUETS

Ø 3.2mm; with Ratchet



**SKU:** 141-030-32

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Catchers, Tourniquets</a>
<b>Aesculap</b>	FB651
<b>KLS Martin</b>	24-556-33
<b>AUTHORS</b>	<a href="#">RUMEL</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Catchers, Tourniquets</a>
<b>Aesculap</b>	FB651
<b>KLS Martin</b>	24-556-33
<b>AUTHORS</b>	<a href="#">RUMEL</a>