

SANTULLI BULLDOG CLAMP

angled; atraumatic

7.0cm (2.8")


SKU: 026-195-07

Categories: [Open Surgery](#), [Hemostats](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| KLS Martin | 24-295-07 |
| AUTHORS | SANTULLI |
| SPECIFICATIONS | angled , atraumatic |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| KLS Martin | 24-295-07 |
| AUTHORS | SANTULLI |
| SPECIFICATIONS | angled , atraumatic |