

## SAROT HAEMOSTATIC FORCEPS

curved; serrated

24.0cm (9.6")




**SKU:** 040-451-24

**Categories:** [Open Surgery](#), [Hemostats](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>KLS Martin</b>	13-433-24
<b>AUTHORS</b>	<a href="#">SAROT</a>
<b>Medicon</b>	16.15.74
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Bronchus Forceps</a>
<b>KLS Martin</b>	13-433-24
<b>AUTHORS</b>	<a href="#">SAROT</a>
<b>Medicon</b>	16.15.74

## SPECIFICATIONS

[curved](#), [serrated](#)