

SAROT HAEMOSTATIC FORCEPS

straight; serrated

24.0cm (9.6")



SKU: 040-450-24

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Bronchus Forceps
Aesculap	BH224
KLS Martin	13-432-24
AUTHORS	SAROT
Medicon	16.14.74
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Bronchus Forceps
Aesculap	BH224
KLS Martin	13-432-24

AUTHORS

[SAROT](#)

Medicon

16.14.74

SPECIFICATIONS

[serrated](#), [straight](#)