

SCHEINMANN EXCHANGEABLE TIP

through-cutting

for BRÜNINGS interchangeable tubes



SKU: 131-485-00

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP226
KLS Martin	39-613-01
AUTHORS	SCHEINMANN
Medicon	73.43.37



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP226
KLS Martin	39-613-01
AUTHORS	SCHEINMANN

Medicon

73.43.37