

SCHEINMANN EXCHANGEABLE TIP

spoon

for BRÜNING'S interchangeable tubes

SKU: 131-484-00

Categories: [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP212
KLS Martin	39-611-01
AUTHORS	SCHEINMANN
Medicon	73.43.27



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Universal handles and tips
Aesculap	OP212
KLS Martin	39-611-01
AUTHORS	SCHEINMANN

Medicon

73.43.27