

SCHMIEDEN LIGATURE CONDUCTOR


curved
 19.0cm (7.6")



SKU: 056-507-19
 Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM906
KLS Martin	20-872-19
AUTHORS	SCHMIEDEN
Medicon	13.32.18
SPECIFICATIONS	curved



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM906
KLS Martin	20-872-19

AUTHORS

[SCHMIEDEN](#)

Medicon

13.32.18

SPECIFICATIONS

[curved](#)