

SCHNIDT TONSIL FORCEPS


curved; open ring
 18.5cm (7.4")



SKU: 132-185-19
Categories: [Open Surgery](#), [Tonsils](#), [Tracheotomy](#), [Tongue](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsils Forceps
Aesculap	BH963
KLS Martin	39-449-19
AUTHORS	SCHNIDT
Medicon	17.21.55



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tonsils Forceps
Aesculap	BH963
KLS Martin	39-449-19
AUTHORS	SCHNIDT

Medicon

17.21.55