

## SCHUMACHER EXCHANGEABLE TIP

through-cutting

for BRÜNINGS interchangeable tubes


**SKU:** 131-492-00

**Categories:** [Goitre](#), [Laryngology](#), [Mouth](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP230
<b>KLS Martin</b>	39-639-01
<b>AUTHORS</b>	<a href="#">SCHUMACHER</a>
<b>Medicon</b>	52.68.57



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP230
<b>KLS Martin</b>	39-639-01
<b>AUTHORS</b>	<a href="#">SCHUMACHER</a>

**Medicon**

52.68.57