

SCISSORS; POINTED; Ø 1.6MM; FR. 5; 34CM

semi-rigid; single action; Luer-Lock
with handle

SKU: 583-229-16LL
Categories: [Endoscopy](#), [Gynecology](#), [GUBBINI/Gyn](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Micro Instruments Hysteroscopy
TYPE	"Semi-Rigid" with LuerLock



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Micro Instruments Hysteroscopy
TYPE	"Semi-Rigid" with LuerLock