


**SCISSORS; POINTED; SINGLE ACTION**

Ø 3.0mm; FL:400mm

with irrigation adapter



**SKU:** 630-461-30

**Categories:** [Endoscopy](#), [Urology](#)

**ZUSÄTZLICHE INFORMATIONEN**

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Nephroscopy</a> |
| <b>TYPE</b>  | <a href="#">Scissors</a>    |

**ZUSÄTZLICHE INFORMATIONEN**

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Nephroscopy</a> |
| <b>TYPE</b>  | <a href="#">Scissors</a>    |