

SEALING CAP BLUE

for trocar sleeve Ø 11.0mm

10 pcs./box; Price per box



SKU: 626-170-90

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------|
| UNIT | PU |
| RANGE | Trocar Systems |
| TYPE | Sheaths |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|--------------------------------|
| UNIT | PU |
| RANGE | Trocar Systems |
| TYPE | Sheaths |