

## SEALING CAP GREEN

for trocar sleeve Ø 12.5mm

10 pcs./box; Price per box



**SKU:** 626-170-11

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>