

SEALING SET Ø 11.0MM; CONSISTING OF:

sealing cap blue; inner sealing yellow

O-ring red;10 pcs./box; Price per box

SKU: 626-171-90

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths