

## SEALING SET Ø 11.0MM; CONSISTING OF:

sealing cap blue; inner sealing yellow

O-ring red;10 pcs./box; Price per box



**SKU:** 626-171-90

**Categories:** [Endoscopy](#), [Laparoscopy](#)

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>