

SEALING SET Ø 12.5MM; CONSISTING OF:

sealing cap green; inner sealing yellow

O-ring; 10 pcs./box; Price per box

SKU: 626-171-11

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths

ZUSÄTZLICHE INFORMATIONEN

UNIT	PU
RANGE	Trocar Systems
TYPE	Sheaths