

SEALING SET Ø 12.5MM; CONSISTING OF:

sealing cap green; inner sealing yellow

O-ring; 10 pcs./box; Price per box



SKU: 626-171-11

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [PU](#)
RANGE [Trocar Systems](#)
TYPE [Sheaths](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [PU](#)
RANGE [Trocar Systems](#)
TYPE [Sheaths](#)