

## SEALING SET Ø 5.5MM; CONSISTING OF:

sealing cap red; inner sealing grey

O-ring red;10 pcs./box; Price per box

**SKU:** 626-171-50

**Categories:** [Endoscopy](#), [Laparoscopy](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">PU</a>
<b>RANGE</b>	<a href="#">Trocar Systems</a>
<b>TYPE</b>	<a href="#">Sheaths</a>