

## SEIFFERT EXCHANGEABLE TIP

toothed

for BRÜNINGS interchangeable tubes

**SKU:** 131-482-00

**Categories:** [Open Surgery](#), [Goitre](#), [Laryngology](#), [Mouth](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP202
<b>KLS Martin</b>	39-604-01
<b>AUTHORS</b>	<a href="#">SEIFFERT</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Universal handles and tips</a>
<b>Aesculap</b>	OP202
<b>KLS Martin</b>	39-604-01
<b>AUTHORS</b>	<a href="#">SEIFFERT</a>