

SEMKEN DRESSING FORCEPS

straight; serrated


15.0cm (6"); width of jaw 1.5 +/-0.2

SKU: 020-040-15
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD216
KLS Martin	12-172-15
AUTHORS	SEMKEN
Medicon	06.30.15
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD216
KLS Martin	12-172-15

AUTHORS

[SEMKEN](#)

Medicon

06.30.15

SPECIFICATIONS

[serrated](#), [straight](#)