

SEMKEN DRESSING FORCEPS

straight; serrated


12.5cm (5"); width of jaw 1.3 +/-0.2

SKU: 020-040-12
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD215
KLS Martin	12-172-12
AUTHORS	SEMKEN
Medicon	06.30.12
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD215
KLS Martin	12-172-12

AUTHORS

[SEMKEN](#)

Medicon

06.30.12

SPECIFICATIONS

[serrated](#), [straight](#)