

## SENN-MILLER RETRACTOR

3 teeth; sharp; double ended

16.0cm (6.4")

**SKU:** 062-180-16

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	BT008
<b>KLS Martin</b>	15-166-16
<b>AUTHORS</b>	<a href="#">MILLER</a> , <a href="#">SENN</a>
<b>Medicon</b>	20.05.08
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	BT008
<b>KLS Martin</b>	15-166-16

**AUTHORS**

[MILLER](#), [SENN](#)

**Medicon**

20.05.08

**SPECIFICATIONS**

[sharp](#), [toothed](#)