


SIMPLEX TUBING CLAMP

with fixation screw

12mm



SKU: 250-400-12

Categories: [Open Surgery](#), [Various Items](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-410-12
AUTHORS	SIMPLEX

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tubing Clamps
KLS Martin	63-410-12
AUTHORS	SIMPLEX