

## SIMPLEX TUBING CLAMP

with fixation screw

12mm



**SKU:** 250-400-12

**Categories:** [Open Surgery](#), [Various Items](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>KLS Martin</b>	63-410-12
<b>AUTHORS</b>	<a href="#">SIMPLEX</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Tubing Clamps</a>
<b>KLS Martin</b>	63-410-12
<b>AUTHORS</b>	<a href="#">SIMPLEX</a>