

## SKILLERN SPHENOIDALIS PROBE

Ø 2.0mm; 16.5cm

**SKU:** 086-702-20

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Probes</a>
<b>AUTHORS</b>	<a href="#">SKILLERN</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Probes</a>
<b>AUTHORS</b>	<a href="#">SKILLERN</a>