

SPENCER-WELLS HAEMOSTATIC FORCEPS

curved; serrated

13.0cm (5.2")




SKU: 040-291-13

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH333 |
| KLS Martin | 13-249-13 |
| AUTHORS | SPENCER, WELLS |
| Medicon | 15.77.13 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH333 |
| KLS Martin | 13-249-13 |

AUTHORS

[SPENCER](#), [WELLS](#)

Medicon

15.77.13

SPECIFICATIONS

[curved](#), [serrated](#)