

SPLINTER FORCEPS

straight; serrated; without Pin

12.5cm (5")

SKU: 020-662-13
 Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-604-12
Medicon	07.56.53
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-604-12
Medicon	07.56.53
SPECIFICATIONS	serrated , straight