

SPLINTER FORCEPS

straight; serrated; with Pin

9.0cm (3.6")

SKU: 020-666-09
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD301
KLS Martin	12-604-09
Medicon	07.80.09
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD301
KLS Martin	12-604-09
Medicon	07.80.09

SPECIFICATIONS

[serrated](#), [straight](#)