

**SPLINTER FORCEPS**

straight; serrated; with Pin  
 10.5cm (4.2")

**SKU:** 020-666-10  
**Categories:** [Open Surgery](#), [Forceps](#)



**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                     |
|----------------|-----------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                                |
| TYPE           | <a href="#">Splinter Forceps</a>                    |
| Aesculap       | BD302                                               |
| KLS Martin     | 12-604-10                                           |
| Medicon        | 07.80.10                                            |
| SPECIFICATIONS | <a href="#">serrated</a> , <a href="#">straight</a> |



**ZUSÄTZLICHE INFORMATIONEN**

|            |                                  |
|------------|----------------------------------|
| UNIT       | <a href="#">pcs.</a>             |
| TYPE       | <a href="#">Splinter Forceps</a> |
| Aesculap   | BD302                            |
| KLS Martin | 12-604-10                        |
| Medicon    | 07.80.10                         |

## SPECIFICATIONS

[serrated](#), [straight](#)