

SPLINTER FORCEPS

straight; serrated; with Pin
 11.5cm (4.6")

SKU: 020-666-11
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD303
KLS Martin	12-604-11
Medicon	07.80.11
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD303
KLS Martin	12-604-11
Medicon	07.80.11

SPECIFICATIONS

[serrated](#), [straight](#)