

SPLINTER FORCEPS

straight; serrated; with Pin
 12.5cm (5")

SKU: 020-666-13
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD304
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD304
SPECIFICATIONS	serrated , straight