

SPLINTER FORCEPS

straight; smooth; without Pin

10.5cm (4.2")

SKU: 020-660-10
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-600-10
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-600-10
SPECIFICATIONS	straight