

SPLINTER FORCEPS

straight; serrated; with Pin
 14.5cm (5.8")

SKU: 020-666-14
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD305
KLS Martin	12-120-14
Medicon	06.02.14
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD305
KLS Martin	12-120-14
Medicon	06.02.14

SPECIFICATIONS

[serrated](#), [straight](#)