

SPLINTER FORCEPS

curved; serrated; with Pin

8.0cm (3.2")

SKU: 020-667-08
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Splinter Forceps](#)
SPECIFICATIONS [curved](#), [serrated](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Splinter Forceps](#)
SPECIFICATIONS [curved](#), [serrated](#)