

SPLINTER FORCEPS

curved; serrated; with Pin

9.0cm (3.6")

SKU: 020-667-09
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD311
KLS Martin	12-605-09
Medicon	07.81.09
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD311
KLS Martin	12-605-09
Medicon	07.81.09

SPECIFICATIONS

[curved](#), [serrated](#)