

SPLINTER FORCEPS

curved; serrated; with Pin

13.0cm (5.2")

SKU: 020-667-13
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD314
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
Aesculap	BD314
SPECIFICATIONS	curved , serrated