

SPLINTER FORCEPS

straight; smooth; without Pin

11.5cm (4.6")

SKU: 020-660-11
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-600-11
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Splinter Forceps
KLS Martin	12-600-11
SPECIFICATIONS	straight