

SPRAY-PROBE

round metal tip type Spray, w. mandrin

1.8mm, length 250cm; single use



SKU: PK-8MREW1802250Q-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)