

STANDARD DRESSING FORCEPS

straight; serrated; narrow

11.5cm (4.6")

SKU: 020-004-11
Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD023
KLS Martin	12-120-11
AUTHORS	STANDARD
Medicon	06.01.11
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Dressing Forceps , Tissue Forceps
Aesculap	BD023
KLS Martin	12-120-11

AUTHORS

[STANDARD](#)

Medicon

06.01.11

SPECIFICATIONS

[serrated](#), [straight](#)