

STANDARD DRESSING FORCEPS

straight; serrated; narrow

18.0cm (7.2")

SKU: 020-004-18

Categories: [Open Surgery](#), [Forceps](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD030 |
| KLS Martin | 12-120-18 |
| AUTHORS | STANDARD |
| Medicon | 06.01.18 |
| SPECIFICATIONS | serrated , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Dressing Forceps , Tissue Forceps |
| Aesculap | BD030 |
| KLS Martin | 12-120-18 |

AUTHORS

[STANDARD](#)

Medicon

06.01.18

SPECIFICATIONS

[serrated](#), [straight](#)